Document ID: 32140582

DX-D 400

Includes the following attached document:

AGFA – DX-D 400 Training and Instruction Protocol (section 9)

Purpose of this Document

This document describes how to:

- Prepare the DX-D 400 Service Application tasks
- Configure the DX-D 400 System
- Train the customer on the DX-D 400 System

Document History

Edition. Revision	Release Date	Changes compared to previous version 1.4
1.5	09-2016	Added information about DR 10s and DR 14s detectors.

Referenced Documents

Document	Reference		
Documents are referenced in the corresponding sections.			

DOCUMENT CONTROL NOTE:
The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.



Manufacturer

Agfa HealthCare N.V.

Published by

Agfa-Gevaert HealthCare GmbH Tegernseer Landstraße 161 D - 81539 München Germany

Copyright © 2016 Agfa HealthCare N.V.

All rights reserved.

Technical modifications reserved.

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert N.V., Belgium, or its affiliates.

All other trademarks mentioned in this document are held by Agfa HealthCare N.V. or the respective owners and are used in an editorial fashion with no intention of infringement.

Nothing contained in this legal notice nor in any text in this document shall be construed as granting by implication, estoppel or otherwise, any license or right to use any of the trademarks, service marks, trade names or logos appearing in this document without the express prior written consent of their respective owner.



WARNING:

Improper operation or service activities may cause damage or injuries.



INSTRUCTION:

- (1) Read the "Generic Safety Directions" prior to attempting any operation, repair or maintenance task on the equipment. Refer to Document ID 11849633.
- (2) Strictly observe all safety directions within the "Generic Safety Directions" and on the product.

LIST OF CONTENTS

1	DX-D 4	100 PRODUCT DESCRIPTION	5
	1.1	Solution Definition	6
	1.2	DX-D 400 Configurations	7
	1.2.1	Wall Stand	12
	1.2.2	Table	12
	1.2.3	Floor mounted Tubestand	13
	1.2.4	Collimators	14
	1.2.5	X-Ray Tube	14
	1.2.6	Generator	15
	1.2.7	Tube / Generator compatibility	15
	1.2.8	Grid	16
	1.2.9	Control Console	18
	1.2.10	Detector	20
	1.2.11	CR Unit	23
2	SAFET	Y DIRECTIONS	24
3	TASK	OVERVIEW	25
	3.1	DX-D 400 analog manual (Config.1 and 2)	25
	3.2	DX-D 400 analog manual (Config.3)	25
4	PREPA	ARATION ACTIVITIES	27
5	SYSTE	EM TAKEOVER FROM FSE	28
6	CONFI	GURATION	29
	6.1	Hard Console – Anatomical Programmer	
	6.1.1	Configuration of the <anatomical programmer=""> (APR)</anatomical>	
	6.1.2	Backup of the Anatomical Programmer File on the Touch Screen Console	
	6.2	Configuring the NX Workstation - DR specific Configuration for DX-D 400 digital manual (D) (Config.3)	31
	6.2.1	Overview of configuration steps	31
	6.2.2	Prerequisites	32

DOCUMENT CONTROL NOTE:
The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.

	6.2.3	Prepare the NX configuration with the NX Offline Config Tool	32
	6.2.4	The default Exam Tree	37
	6.2.5	Importing the Exam Tree - local default or customized	37
	6.2.6	DX-D 400 specific Exam Tree configuration	38
	6.2.7	Compose exam	44
	6.2.8	Arranging XRG parameters / User Interface configuration	45
	6.2.9	Enable/disable Create DR Sequence button	46
	6.2.10	Activate the configuration	46
7	CUSTO	OMER TRAINING AND HANDS-ON ASSISTANCE	47
	7.1	Training Prerequisites and Preparation	47
	7.2	Training Checklist – General Topics	48
	7.3	Training Checklist – X-Ray Unit	49
	7.3.1	Training Checklist – Manual (Touch Button) Console	51
	7.3.2	Training Checklist – Soft Console	53
	7.4	Detector	54
	7.4.1	Wireless Detector - DX-D 30C/35C	54
	7.4.2	Wireless Detector - DX-D 40C/G and DX-D 45C/G	55
	7.4.3	Wireless Detector DR 14C/G and DR10C	56
	7.4.4	Checklist for Customer Training "DR Detector Calibration"	58
	7.5	Checklist for Customer Training "NX Workstation"	59
8	DX-D 4	100 GLOSSARY	63
9	AGFA -	– DX-D 400 - TRAINING AND INSTRUCTION PROTOCOL	64

1 DX-D 400 Product Description

The DX-D 400 floor mounted X-ray System is an X-ray imaging system used in hospitals, clinics and medical practices by radiographers, radiologists and physicists to make, process and view static X-ray radiographic images of the skeleton (including skull, spinal column and extremities), chest, abdomen and other body parts on adult, pediatric or neonatal patients.

Applications can be performed with the patient in the sitting, standing or horizontal position.

The DX-D 400 X-ray unit can be used in the General Radiology:

- Hospitals / Hospital Satellites
- Imaging Centers
- Emergency work in smaller hospitals / imaging centers
- Orthopedic Clinics/Practices
- Private Practices



Figure 1

Limitations:

The DX-D 400 system is not intended for Mammography applications.

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.

1.1 Solution Definition

Depending on the installed configuration the DX-D 400 system can be used either for conventional film/screen (Config.1 or 2), for CR (Config.1 or 2) or as an integrated digital X-ray system (Config.3).

It consists of the following components:

- Floor mounted Tube stand
- · Wall Bucky Stand with cassette tray or fixed detector
- Table with cassette tray or fixed detector
- X-ray tube with tube head display and collimator
- High frequency generator with manual (touch button) console or soft console (touch screen)
- NX workstation depending on configuration




Figure 2

DOCUMENT CONTROL NOTE:

Solution Options:

- Wet film
- CR unit (DX-G, DX-M, CR10-X, CR12-X or CR15-X / CR30-X or CR-30)
- Hardcopy printer

DX-D 400 System Accessories:

- Compression band
- Hand grips for table
- Hand grips for wallstand
- Hand support for wallstand
- Lateral cassette holder 24 x 30 cm
- Lateral cassette holder 35 x 43 cm
- DX-D 400 VACUDAP 2004 (External Dose Area Meter)
- DX-D 400 VACUDAP 2004 (External Dose Area Meter) integrated

1.2 DX-D 400 Configurations

The DX-D 400 is available in three different basic configurations. The components depend on the configuration:

- as a non-integrated analog modality with a hard console without software integration and for CR / Film only:
 - o DX-D 400 analog manual, (Configuration 1)
 - DX-D 400 analog manual (D), (Configuration 2)

and

- as an integrated Digital Radiography X-ray system with software integration and detector(s):
 - o DX-D 400 digital manual (D), (Configuration 3)

Details see table on next page.

Configuration*	Table	Bucky in the wall stand	Tubehead Display	Collimator	Control	Image Processing
DX-D 400 analog manual		CR	Analog	Manual or	Hard console	Wet film, NX workstation
DX-D 400 analog manual (D)	Fixed or elevating	r Grand	automatic	with touch screen	or other CR workstation	
DX-D 400 digital manual (D)		DR	Digital	Manual	Soft console	NX workstation

* Meaning of the configuration names:

- analog: CR - Computed radiography (requires digitizer with cassettes)

- digital: DR – Digital radiography (has two fixed detectors built in)

- manual: Not motorized

- (D): Digital tubehead display

(Config.1)




Figure 3

- Floor mounted tubestand with analog tubehead display
- Wallstand with cassette bucky
- non tilting or tilting
- left or right loading
- optional with Automatic Cassette Size Sensing
 - Floating table
- non elevating or elevating
- with cassette tray
- optional with Automatic Cassette Size Sensing
 - Manual or automatic collimator
 - Hard Console
- optional with APR (Automatic Programmer)
- optional with touch screen incl. APR

(Config.2)




Figure 4

- Floor mounted tubestand with digital tubehead display for angles and SID
- Wallstand with CR bucky
- Non tilting or tilting
- left or right loading
- optional with Automatic Cassette Size Sensing
 - Floating table
- non elevating or elevating
- with cassette tray
- optional with Automatic Cassette Size Sensing
 - Manual or automatic collimator
 - Hard Console
- optional with APR (Automatic Programmer)
- optional with touch screen incl. APR

(Config.3)




Figure 5

- Floor mounted tubestand with digital tubehead display for angles and SID
- Wallstand with DR bucky
 - non tilting or tilting
- Floating table with DR bucky
 - non elevating or elevating
- Manual collimator
- NX with touch screen monitor
- Agfa Softconsole
- DX-D Retrofit Box/DR Generator Sync Box and up to 3 detectors (different combinations)

1.2.1 Wall Stand

Available as non tilting and manual tilting wall stand:

- Vertical bucky assembly
- All standard cassette sizes from 18 to 43 cm are accepted for CR
- Chin rest support at the upper side
- Panel shows AEC detector areas and cassette film size
- Loading of Cassette in Wall Bucky Stand can be adapted at right or left side as per installation requirements.
- Wallstand with fixed detector can be adapted for movement at right or left side.



Figure 6

1.2.2 **Table**

- Maximum patient weight: 350 kg (771 lbs)
- Floating table top
- Non elevating or elevating
- For standard cassette sizes from 18 to 43 cm (Config.1 and 2)



Figure 7

1.2.3 Floor mounted Tubestand

- Support for the Tube-Collimator Assembly with Column Control Panel for:
 - Horizontal movement of column
 - o Vertical movement of the Tube-Collimator Assembly
 - o Rotational Movement of the Tube-Collimator Assembly (0°, 90° and -90°)
 - Transverse movement of the Tube-Collimator Assembly

Config.1 • Analog tubehead display



Figure 8

Config.2-3 • Digital tube head display

 Position data of the Tube-Collimator Assembly (Angle and SID)



Figure 9

1.2.4 Collimators

Config.1-3 • Manual collimator:

- Collimator controls consist of a push-button to switch on the Collimator Lamp
- Two knobs to open or close the internal blades of the Collimator.



Figure 10

Config.1-2 • Automatic Collimator ACSS (automatic cassette size sensing) for CR only

- A digital function display.
- Memorization of the two formats that are more frequently used.
- 5 adjustable fixed fields to be used with one detector (flat panel).
- ACSS Automatic Cassette Size Sensing (for CR only)
- Additional filtration (pushbutton selection)
 - 1mm Al+0.1mm Cu,
 - 1mm Al+0.2mm Cu
 - 2mm Al (1mm Al optional)



Figure 11

1.2.5 X-Ray Tube

X-Ray Tube	Anode angle	Tube Focus / Power	Heat Capacity	
			(kilo Heat Units)	
Toshiba E7884X (default)		0.6 mm / 20 KW 1.2 mm / 50 KW	300 kHU	
Toshiba E7252X	12°	0.6 mm / 27 KW 1.2 mm / 75 KW	300 kHU	
Toshiba E7254FX	12	0.6 mm / 40 KW 1.2 mm / 102 KW	400 kHU	
Toshiba E7869X		0.6 mm / 40 KW 1.2 mm / 100 KW	600 kHU	

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.

1.2.6 **Generator**

Configurations with 50, 64 and 80 kW generators are available.

1.2.7 **Tube / Generator compatibility**

50kw Generator (Low and High Speed)	64kW Generator	80kW Generator
	E7252X,	E7254FX
All	E7254FX	E7869X
	E7869X	E7009X

1.2.8 **Grid**

A variation of grids is available. Grids fixed in the bucky are always moving grids, removable grids are always static (not motorized). The selection of grids depends on the different DX-D 400 versions and the configuration.

Wallstand:

- 103 LINES, 12:1, FFD 1.5 M (CR, film) → non removable grid, moving grid
- 103 LINES, 12:1, FFD 1 M (CR, film) → non removable grid, moving grid
- 103 LINES, 12:1, FFD 1.8 M
- → non removable grid, moving grid
- 132 LINES 10:1 FFD 1M
- 132 LINES 10:1 FFD 1.5M
- 132 LINES 10:1 FFD 1.8M
- 215 LINES, 12:1, FFD 1 M
- 215 LINES, 8:1, FFD 1.3 M
- 215 LINES, 12:1, FFD 1.5M, delivered without frame
- 215 LINES, 12:1, FFD 1.8M, delivered without frame

Table:

- 103 LINES 12:1 FFD 1M
- → non removable grid, moving grid
- 130 LINES 8:1 FFD 1M
- 215 lines, 12:1, FFD 1 M
- 132 lines, 10:1, FFD 1 M

Recommended Grid – Detector combination

	DX-D 10	Varian Fixed	DX-D 30 DX-D 35	DX-D 40 DX-D 45	DR 14s DR 10s
215 LINES 8:1 FFD 1.3M CARBON COVER	✓	✓			
215 LINES 12:1 FFD 1.5M CARBON COVER	✓	✓			
215 LINES 12:1 FFD 1.8M CARBON COVER	√	✓			
215 LINES 12:1 FFD 1 M CARBON COVER	√	✓			
132 LINES 10:1 FFD 1M CARBON COVER			✓	✓	
132 LINES 10:1 FFD 1.5M CARBON COVER			✓	✓	
132 LINES 10:1 FFD 1.8M CARBON COVER			✓	✓	

Refer to the Agfa website for specifications on the anti-scatter grids that have been found compatible with the system and the DR Detectors.

http://www.agfahealthcare.com/global/en/library/overview.jsp?ID=54332498

1.2.9 Control Console

Manual (touch button) Console

DX-D 400 All controls, indicators and displays are functionally (*Config.1*) grouped.

DX-D 400 (Config. 2)

- Power on/off
- Workstation selection
- AEC Module
- Focal spot
- Exposure indicators
- Radiographic values
- Anatomical Programmer
- Exposure technique:
 - kVp with AEC operations one point technique
 - o kVp and mAs two point technique
 - kVp, mA and Exposure Time three point technique



Figure 12



Figure 13

- Anatomical Programmer (APR):
 - With factory pre-programmed anatomical views for automatic selection.
 - For six patient sizes (three adults and three children)
 - The operator may manually modify all the original APR techniques and store them for later use.



Figure 14

X-Ray Soft Console

DX-D 400 • (Config. 3)

- Located on the right-hand side of the NX wide screen Monitor (See red rectangle below.)
- Provides controls for e.g. position information, SID, angle and rotation of tube
- Generator Controls for adjustment of exposure parameters
- Operational Status Bar to show the device status, DAP display and error messages






Figure 15

Figure 16: Generator Control

1.2.10 **Detector**

With release of this document the following detectors are supported:

Integrated detector	
 Fixed in table and/or wallstand ~ 43 x 43 cm 	
DX-D 10	
Portable	
Cassette sized	1 00000
• ~ 35 x 43 cm	
DX-D 30C	
Wireless	1
Cassette sized	DX-D300
• ~35 x 43 cm	
DX-D 35C	att
Wireless	1 - 1
Cassette sized	
• ~ 31 x 38 cm	
DX-D 40C/G	
Wireless	
Cassette sized	0x-0 40C
• ~ 37 x 44 cm	
DX-D 45C/G	
Wireless	
Cassette sized	acces /
• ~ 29 x 35 cm	
DR 14s C/G	
Wireless	MUSICA
Cassette sized	0104
• ~35 x 43 cm	OR 141
DR 10s C	
Wireless	- MOUNTER
Cassette sized	1 + 1
• ~ 24 x 30 cm	on son of the son of t

Figure 17

1.2.10.1 **Basic Detector Data**

	Fixed Detectors	DX-D 10	DX-D 30C	DX-D 35C	DX-D	40C/G	DX-D 45C/G
		01000	DX-030C		2. 040	100	4000
Receptor Type	Amorphous Silicon						
Scintillator	Csl or GOS	Csl or GOS	Csl	Csl	Csl or	GOS	Csl or GOS
Data transmission	Cable	Cable	Wireless	Wireless	Wire	less	Wireless
Resolution	3.6 lp/mm (CsI) 2.8 lp/mm (GOS)	3.6 lp/mm (CsI) 2.8 lp/mm (GOS)	4 lp/mm	4 lp/mm	3.1 lp/mm (CsI) 2.8 lp/mm (GOS)		
Energy Range - Standard				40 - 150 kV			
Dynamic Dose Range	Csl: 0.3 microgra GOS:0.3 microgra		0.03- 50 ו	microgray	Max. Exposure Level 90 microgray		
Size	42.7 x 42.7 cm (16.8 x 16.8 inch)	42.7 x 35.6 cm (16.8 x 14.0 inch)	46 x 38.4 cm (18.1 x 15.1 inch)	30,7 x 38,4 cm (12.1 x 15.1 inch)	36,8 x ² (14,49 x 1		approx. 28,7 x 35,0 cm (11,3 x 13,8 inch)
Effective Image Area	42.4 x 42.4 cm (16.7 x 16.7 inch)	42.4 x 35.3 cm (16.7 x 13.9 inch)	35.0 x 42.6 cm (13.8 x 16.8 inch)	27,4 x 35 cm (10,9 x 13,8 inch)			25,2 mm x 31,6 mm (9,92 x 12,45 inch)
Weight		approx. 3.9 kg (8.6 lbs.)	3.4 kg (7.5 lbs) (incl. battery)	2.3 kg (5.1 lbs) (incl. battery)	DX-D 40G 3,3 kg (7.3 lbs) (incl. battery)	DX-D 40C 3,4 kg (7.5 lbs) (incl. battery)	DX-D 45C/G 2.2 kg (4,9 lbs) (incl. battery)

DOCUMENT CONTROL NOTE:
The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.

	DR 14s (Pixium 3543 EZ)	DR 10s (Pixium 2430 EZ)			
	DE SEL	on the St. St.			
Receptor Type	Amorphous Silie	con			
Scintillator	Csl or GOS	Csl			
Data transmission	Wireless				
Resolution	3,4 lp/mm	3,4 lp/mm			
Energy Range - Standard	40 - 150 kV				
Dynamic Dose Range	CsI : max linear dose ≥ 50μGy;, clipping dose > 85 μGy GOS: max linear dose ≥ 75μGy; clipping dose > 125 μGy	max linear dose ≥ 50μGy; clipping dose > 85 μGy			
Size	approx. 38.5 x 46 x 1.6 cm (15 x 18 x 6.3 inch)	approx. 26.9 x 32.9 x 1.6 cm (10.6 x 12.6 x 6.3 inch)			
Effective Image Area	33.7 cm x 41.4 cm (13.3 x 16.3 inch)	22.2 cm x 28.4 cm (8.7 x 11.2 inch)			
Weight	DR 14s C: 3.0 kg (6.6 lbs) DR 14s G: 2.8 kg (6.2 lbs)	1.6 kg (3.5) lbs (incl. battery)			

DOCUMENT CONTROL NOTE:
The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.

1.2.10.2 Possible detector combinations

- Two fixed detectors in radiographic table and wall stand OR
- Fixed detector in the wallstand and cassette size detector in the table OR
- Single cassette size detector to be swapped between table and wallstand

1.2.11 CR Unit

Supported Agfa Digitizers:

• CR10-X

DX-G

• CR12-X

• DX-M

• CR15-X / CR30-X

References for Digitizer training preparation in the Agfa HealthCare Library:

Document Type	Agfa HealthCare Library	
CR 10-X, CR 12-X, CR 15-X User Manual 2491		
CR 30-X Changing from analog to digital UM 2397	Computed Radiography → CR Digitizers → User Manual	
CR 30-X/CR 30-Xm User Manual 2386		
DX-G/DX-M User Manual 2321		

2 Safety Directions



WARNING:

Improper operation or service activities may cause damage or injuries.



INSTRUCTION:

- (1) Read the "Generic Safety Directions" prior to attempting any operation, repair or maintenance task on the equipment. Refer to Document ID 11849633.
- (2) Strictly observe all safety directions within the "Generic Safety Directions" and on the product.

Please refer to product-specific safety directions available in the different User Manuals applicable for the DX-D 400 system.

User Manuals are part of the delivery.

They can also be downloaded from the Agfa HealthCare Library.

Available Documentation	Path
DX-D 400 User Manual 0230	Direct Radiography →
DX-D 400 X-Ray System User Manual 0232	DX-D 400 → User Manual
DX-D Software Console User Manual 0189	Direct Radiography → DR Components → User Manual
DX-D 30C, DX-D 35C User Manual 0197	Direct Radiography →
DX-D 40C, DX-D 40G User Manual 0290	DR Detectors → User Manual
DX-D 45C, DX-D 45G User Manual 0292	
DX-D 10G, DX-D 10C, DX-D 20G, DX-D 20C User Manual 0129	
DR 10s Detector User Manual 0351	
DR 14s Detector User Manual 0350	
DX-D DR Detector Calibration Key User Manual 0134	

3 Task Overview

The following tables give an overview of basic Service Application tasks.

3.1 DX-D 400 analog manual (Config.1 and 2)

Overview

Details see section	Activity	Description		Required Time
4	Preparation Activities	Prepare Service Application Tasks after Customer Contract Signature (organize visit).		2 hours
5	System Takeover	System Takeover of the system from Field Service Engineer after installation		15 min
6	Configuration Tasks (= system setup)	Perform default configuration		2 hours
7	Customer Training and Hands-on Assistance on system level	Radiographer training on X-ray system	Operation of X-ray system	3 hours

3.2 DX-D 400 analog manual (Config.3)

Overview

Details see section	Activity	Description		Required Time
4	Preparation Activities	Prepare Service Application Tasks after Customer Contract Signature (organize visit).		2 hours
5	System Takeover	System Takeover of the system from Field Service Engineer after installation		15 min
6	Configuration Tasks (= system setup)	Perform default configuration		2 hours
7	Customer Training and Hands-on Assistance on system level		Operation of X-ray modality	
	on system level	Radiographer training on NX,	Operation and workflow of NX workstation	0.1
		X-ray modality and integrated detector	System workflow	2 hours
			Calibration of detector	
		Assist customer in system operation and answer questions to finalize tasks		4 hours

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.

Prerequisites

Prerequisites to keep the timing for the basic application:

The Sales representative has to inform the customer about the activities provided by Application as listed in the table above. No customization is foreseen. Extra activities need to be ordered separately by the customer.

Restrictions DX-D 400 analog manual (Config.1) and DX-D 400 analog manual (D) (Config.2):

- No adaptation of the Anatomical Programmer (APR) with exam tree and X-ray parameters.
- No hands-on session with patients is foreseen.

DX-D 400 digital manual (Config.3):

- 1. No adaptation of the DX-D 400 local exam tree and X-ray parameters
- 2. No protocol codes configured
- 3. No image quality adaptations including taste settings Musica²

Preparation Activities 4

DX-D 400 - Checklist for Preparation Activities			
Site:	Location:	Date:	
Goal:			
The following list sun Service Application to	nmarizes the activities to prepare the DX asks.	-D 400	
Topics:		Done	Reference
Preparation Activities	es		
Identify persons and which are involved in	representatives of the different groups the project.		Not applicable
together with the cust	/ customization tasks for the solution tomer. Keep in mind the timing for the tra activities need to be ordered stomer.		
Observe if extra appli have been ordered.	ication relevant tasks are necessary and		
(APR) strongly depen	configure the anatomical programmer ands on its size. Estimate the additional of the programs, if the customer does t.		
Tree, additional time considered. This addicustomer. Note down requested	not accept the Default DX-D 400 Exam for exam tree configuration has to be itional effort should be paid by the I adaptations, e.g. add exposures, renaming of exposures.		
Prepare the exam tre	e with the NX Offline Config Tool.		
Analyze the time plan	according to application relevant tasks.		
	nd guidelines on regional, country and ve to be filtered with respect to the ect.		
Analyze the Custome	er Training Requirements.		
Define together with t training should take p	the customer when and where the place.		
Hand over the filled-in Project Manager.	n Training and Instruction Protocol to the		
End			

DOCUMENT CONTROL NOTE:
The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.

System Takeover from FSE 5

DX-D 400 - Checklist for System Takeover			
Site:	Location:	Date:	
Goal:			
	opics to be checked with the Field -up of the system is complete.	d Service	Engineer to
Topics:		Done	Reference
System Takeover Activities			
Check that all components or present and technically set-up			Not applicable
Check that the DX-D 400 Sys on site (e.g.: on Documentation delivery).	stem User Manuals are available on DVD, which is part of the		
Check if a backup has been of files:	created for the APR configuration		
APRBckUp.ini			
APR.ini			
Check also for the location ar	nd the naming of these files:		
If no backup has been perfori	med, refer to section 6.		
Check which kind of detector Cesium lodide) is installed	(Gadolinium Oxysulfide or		
Have the sensitivity settings to ones?	een communicated? Which		
Medium: μGy			
Low and High sensitivity dose calculated.	e values are mathematically		
Get informed about any open unspecified settings, values of	issue: e.g. missing components, or results.		
Make sure that any open tech escalated according to the Ag			
End			

6 Configuration

6.1 Hard Console – Anatomical Programmer

6.1.1 Configuration of the <Anatomical Programmer> (APR)

On the touch console an exam tree can be configured in the <Anatomical</pre>
Programmer> (APR).

The available default configuration can be changed either by the customer or service.

For detailed information of configuration possibilities refer to: 0232 DX-D 400 X-Ray System User Manual.

6.1.2 Backup of the Anatomical Programmer File on the Touch Screen Console

It is strongly recommended to make a backup of the APR settings after installation and configuration of the system.

Two different files are available on flash card in the touch console:

- APRBckUp.ini → APR back up file
- APR.ini → APR working file



NOTE:

Make sure a card reader is available.

- (1) Switch off the system.
- (2) Unscrew and remove the generator cable.



Figure 18

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.

(3) Unscrew and remove rear cover (4 screws).



Figure 19

(4) Take the flash card out.

File location on flash card is:

C:\Nuevo\program files\TPC\APR



Figure 20

- (5) Make a copy of the files:
 - APRBckUp.ini
 - o APR.ini
- (6) Rename the copied files to, e.g.:
 - o defaultAPRBckUp.ini
 - o defaultAPR.ini
- (7) Return the flash card in the console. Mind its correct position.
- (8) Mount the rear cover and connect the generator cable.
- (9) Switch the system on.

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.

6.2 Configuring the NX Workstation - DR specific Configuration for DX-D 400 digital manual (D) (Config.3)



Required time for a DX-D specific NX configuration:

Approximately 30 minutes if a local and customized Exam Tree is used. Approximately 4 hours and more have to be estimated only for Exam Tree configuration and customization if the standard default Exam Tree is used.

6.2.1 Overview of configuration steps

DX-D 400 - Checklist for NX Configuration			
Site: Location:	Date:		
Task		Section	
Prepare the NX configuration with the NX Offline Config Tool		6.2.3	
Configure the X-Ray Device		6.2.3.1	
Configure the DR Detector		6.2.3.2	
Configure the DR recovery procedure		6.2.3.3	
Exam Tree Configuration		6.2.4	
Import the local default exam tree for DX-D 400 or the customized exam tree.		6.2.5	
Perform the DX-D 400 specific configuration of the Exam Tree.		6.2.6	
Configure the exposure settings.		6.2.6.1	
Configure film orientation (optional).		6.2.6.2	
Configure the modality settings.		6.2.6.3	
Configure Add Image Pane in Compose Exam.		6.2.7	
Arrange XRG Parameters / User Interface Configuration.		6.2.8	
Enable/disable Create DR Sequence button		6.2.9	
Activate settings.		6.2.10	

DOCUMENT CONTROL NOTE:

6.2.2 Prerequisites

- NX Workstation is installed.
- XRDI software is installed.
- CDI software and DR Detector-related software is installed and configured.

evice Configuration

- Default DX-D 400 Exam Tree is available.
- Licenses are enabled.

6.2.3 Prepare the NX configuration with the NX Offline Config Tool

6.2.3.1 Configure the X-Ray Device

- (1) To add an X-ray device:
 - Click the **New** button.
 - Type in a reasonable name.
 - Click **Upload** to select the corresponding X-ray device model file.



Figure 21

(2) Mark the model file.
Click **Open** and then **OK**.




Figure 22

DOCUMENT CONTROL NOTE:

The following table shows a list of available model files:

Available system	X-Ray device model file to be selected
DR 400 Type 5520/100	Spellman.DR400.xml
DR 400 Type 5520/200	Agfa.DR400.xml
DX-D 100 system with a 20 kW generator	Sedecal.DXD100_20kW.xml
DX-D 100 system: All except 20 kW generator	Sedecal.DXD100.xml
DX-D 300 system	Sedecal.XPLUS.xml
DX-D 400, with or without Retrofit Box/DR Generator Sync Box	SedecalAM.xml
DX-D 500 system	Siemens.Polydoros+DR.xml
DX-D 600 system	Sedecal.DXD600.xml
Livingston Mammo generator	Livingston.xml
Planmed Mammo generator	Planmed.Sophie.xml
Quantum solution CR and DR	Quantum.Techvisison.QBus+DR.xml
Quantum solution only CR	Quantum.Techvisison.QBus.xml
Retrofit system with Sedecal generator connection (XRDI connection)	SedecalAM.xml
Retrofit system with Siemens generator connection (XRDI connection)	SiemensAM.xml
Retrofit system without XRDI connection between NX and generator.	DRPanelConnectionOnly.xml
Siemens DXSi (XRDI 10)	GENIT generator.xml
Siemens DXSI Solution (from XRDI 15 and higher)	Siemens.DXSi.xml.xml
Siemens Mammo	Siemens.Mammomat.xml



IMPORTANT:

Do NOT use the following model files:

• RetrofitBox.xml: R&D internal model file. Not relevant for service.

Sedecal.DXD400.B.xml: Not used.
 Sedecal.DXD600.xml: Not used.
 Sedecal.DXD600.B.xml: Not used.

6.2.3.2 Configure the DR detector

- (1) To add a DR detector:
 - Select the X-Ray device.
 - Click Add.



Figure 23

- Highlight the required DR detector model.
- Click Open and click OK.




Figure 24

(2) In Detector Settings use the dropdown menu for the **Serial Number** and select:

Wall, Table or Portable




Figure 25

(3) Type in the **Detector Nickname**, e.g. Wall, Table, Large, Small ...




Figure 26



NOTE:

The DX-D Retrofit Box/DR Generator Sync Box settings do NOT need to be configured in the NX Offline Configuration tool.




Figure 27

6.2.3.3 Configure the DR Recovery Procedure

If an exposure is no longer linked to a thumbnail, it will come in on NX as emergency exposure. The image processing applied to this kind of exposure is defined via Exposure Type.

To configure the DR Recovery Procedure at least one DR exposure must be set up in the exam tree.

- (1) In **Devices** mark X-Ray Device.
- (2) Select Age Group.
- (3) Select an Exam Group.
- (4) Select an Exposure Type.




Figure 28

6.2.4 The default Exam Tree



IMPORTANT:

- The default factory Exam Tree for GenRad that is delivered on each new NX workstation does not include the modality settings, e.g. the exposure settings and receptor type.
- It is the responsibility of each country/region to optimize the configuration of the default DX-D 600 Exam Tree according to local requirements and to re-use this configuration to avoid unnecessary repetition of work performed on site.

6.2.5 Importing the Exam Tree – local default or customized

To import either the local default exam tree or the customized exam tree, follow this procedure:

- (1) Open the Configuration Tool.
- (2) Go to: Exam Tree Configuration
- (3) Select: Add Exams from Import...
- (4) Go to the directory where either the local default exam tree or the customized exam tree is stored, e.g.: C:\AGFA\Healthcare\WX\bin\Factory Settings
- (5) Open the prepared exam tree from that folder.
- (6) Click **OK**.
- (7) Select the exams to be included and click **Import**.



NOTE:

After the import, all exposures are visible in the exam tree, regardless of the device configuration. However, in the NX GUI the user will only see those exposures that are relevant.

Example:

If you import an exam tree with CR exposures, but no digitizer is configured, the CR exposures are visible in the Configuration Tool. In contrast, no CR exposures will be visible in the NX GUI after activation.

DOCUMENT CONTROL NOTE:

6.2.6 DX-D 400 specific Exam Tree configuration



NOTE:

For the procedure of how to configure an exam tree on NX, refer to the Service Manual and Key User Manual:

- Chapter 4, Installation and Configuration Appendices of the NX Service Manual
- Chapter 6, Configuring the Examination Tree of the NX Key User Manual

6.2.6.1 Configuration of the exposure settings

It is possible to configure the Auto Rotate (180 °) function on exposure level to prevent that the customer has to rotate the image manually. This may appear for configurations like right and left exposures for e.g. hand, wrist and elbow.

A target exposure index can also be fixed here, instead of using the dose monitoring tool.

- (1) Go to: Exam Tree Configuration
- (2) Select an exposure.
- (3) Check the Auto Rotate 180 degree function.




Figure 29

- (4) Select the Cassette orientation.
 - For fixed detectors always choose orientation Portrait.
- (5) Fix a Target Exposure Index (optional).

DOCUMENT CONTROL NOTE:

6.2.6.2 Configuration of film orientation

On CR systems the film orientation is determined by the cassette orientation. On DR systems with <u>fixed detector</u> the orientation of the detector is fixed to portrait. Therefore the cassette orientation is not applicable for the film orientation.

In the configuration part **Printer Settings** of the exposures the configuration option **Film Orientation** is offered:

- Automatic:
 - The Region of Interest (ROI) is automatically detected and acts as reference for the film orientation (landscape or portrait is automatically chosen).
- Landscape:
 - The film orientation is always set to landscape.
- Portrait:

The film orientation is always set to portrait.




Figure 30



NOTE FOR OPTION FILM ORIENTATION:

- Because of anatomical variation the selected film orientation does not always deliver the optimal result for the printout.
- Inform the user about this behaviour and the possibility to correct the film orientation in the **Editing** pane of the NX user interface.

6.2.6.3 Configuration of modality settings

- (1) Go to: Exam Tree Configuration
- (2) Select an exposure.
- (3) In the **Modality Settings** pane go to section **General settings**.

The view with ViewName "1" opens. The fields for which a configuration activity is required are marked with a red rim.




Figure 31: Example screenshot

(4) Select the correct modality position depending on system configuration, e.g. DR Table, Wall fixed DR detector, Wall DR portable, Table DR portable, Free DR portable, CR Wall, CR Table, CR Free, ...



Figure 32: Example screenshot

DOCUMENT CONTROL NOTE:

(5) Select the **Detector Nickname**, if more than one portable detector is configured. Link the exposure to the detector which will be mainly used for this exposure.




Figure 33

(6) Enter the kVp.

Enter the mAs

or

AEC sensitivity, chambers, density.

Possible ms for fixed detector, DX-D 10 and DX-D 20: 550 or 1000 ms Possible ms for DX-D 30C/DX-D 35C: 1000 or 3000 ms Possible ms for DX-D 40/DX-D 45: either 1000 or 3000 ms (time is defined during installation)



Figure 34: Example screenshot



NOTE:

Optional and in combination with AEC the mA and max. mAs can be configured.

For more information refer to the Application Note – NX configuration items in DR systems max mAs and mA (Document ID: 43122034).

(7) Select the appropriate Focus, Grid and, if necessary, the PreFilter.




Figure 35

(8) Configure the **Exposed Area (mm)** (automatic Collimation).



Figure 36



NOTE:

The default programs 20, 24 and 28 are relative positions with defined SIDs, but not defined wall stand height. The tube will always stop perpendicularly to the detector in the wall stand.

The default Auto-Position programs can only be adapted via the Service Tool in the software of the Tubehead display.

Procedures are trained in the DX-D 600 System Application Training. For booking of Instructor Based Training refer to the Academy Learning Platform (ALP), Hyperlink: https://healthcare.agfa.net/irj

(9) Repeat steps (2) to (9) for all exposures.



IMPORTANT:

With the tool *Edit Exposure Modality Settings* in the exam tree configuration it is possible to adapt single or various exposures.

See Application Note - NX 2.0.8900 - NX 3.0.8900 Configuration Improvements on Agfa HealthCare Library: Document ID:47340342




Figure 37




Figure 38

DOCUMENT CONTROL NOTE:

6.2.7 Compose exam

- (1) Go to: User Interface
- (2) In Examination Settings Select Configure Add Image Pane.
- (3) Open the tab Compose Exam.
- (4) Place the exposure thumbnails for all exams and age groups at the preferred position.




Figure 39

6.2.8 Arranging XRG parameters / User Interface configuration

- (1) Go to: User Interface Configuration → Examination Settings → Configure XRG Parameters
- (2) Open the **XRG Parameters** screen to set the displayed attributes.




Figure 40: Example screenshot

- (3) Select a value for **Position**.
- (4) Translate the **Label** entry to local language, if necessary.
- (5) In **Content** select the XRG Parameter to be displayed.

6.2.9 Enable/disable Create DR Sequence button

- (1) Go to: General Configuration → Workflow Management → Identification
- (2) Enable/disable the **Create DR sequence button** according to the customer workflow.




Figure 41

6.2.10 Activate the configuration

Apply the configuration settings. In the Configuration Tool menu select:

File → Activate Configuration

7 Customer Training and Hands-on Assistance

7.1 Training Prerequisites and Preparation

For a successful training the following assumptions should be met:

- DX-D 400 is fully installed, configured and functional.
- Date of training is defined.
- List of persons to be trained is available.
- The recommended maximum number per training group is 4 participants.
- No interruptions during training and hands-on assistance may occur.

Training Checklist – General Topics 7.2

	DX-D 400 - Checklist for Customer Training				
Site	: Location:	Date:			
Afte	Goal: After the training the customer should understand the workflow of the DX-D 400 and know how to use the documentation.				
Topi	cs:	Done	Reference		
Gen	eral Topics				
•	Intended Use		0232 DX-D 400 User		
•	Available User Manuals		Manual		
•	System Overview				
•	"Emergency Stop" button(s)				
•	Safety directions				
•	Start-up and Switch off Procedure DX-D 400		0232 DX-D 400 User		
	o Start-up Routine		Manual		
	o Shut-down Routine				
	o Rules for detector warm-up				
	 X-ray tube warm-up procedure 				
•	Cleaning and disinfection				
•	Periodic Maintenance				
End					

Training Checklist – X-Ray Unit 7.3

	DX-D 400 - Checklist for Customer Training			
Site) :	Location:	Date:	
Goa				
		participant knows the different controls of the DX-D 40)O.	I
Тор			Done	Reference
Floo	or Mo	unted Tube Stand		
•	DX- or	-D 400 Analog Control Panel		0232 DX-D 400 User
	_	-D 400 Digital Control Panel		Manual
	0	Movements of column		
	0	Movement of the Tube-Collimator Assembly		
	0	Light indicators		
	0	Angle and SID display (digital control panel only)		
• (Colun	nn Rotation Control		
RAD) Tab	le		
•	Fixe	ed Height Table / Elevating Table		0232 DX-D 400 User
	0	Emergency switch		Manual
	0	Circuit breaker		
	0	Table top lock switch		
	0	Bucky handle (CR)		
	0	Cassette tray (CR)		
•	Fixe	ed Height Table		
	0	Table - top brake pedal (longitudinal and transverse motion)		
•	Ele	vating Table		
	0	Table - top movement pedals (longitudinal and transverse motion, up and down)		
RAD) Wal	I Stand		
•	Bud	cky table top		0232 DX-D 400 User
	0	Grooves for standard accessories		Manual
	0	Chin rest support		
	0	AEC detector areas		

		DX-D 400 - Checklist for Customer T	raining	9	
Тор	Topics: Done Reference				
RAI) Wa	II Stand (continued)			
•	Co	ntrols for			
	0	Vertical lock of bucky assembly			
	0	Bucky assembly rotation			
	0	Manual tilting of bucky			
•	Cas	ssette tray (CR)			
Gric	i		•		
•	Kin	d of grid used in table and wallstand		0232 DX-D	
•	Ge	neral rules for using grids, e.g. SID, Range of SID		400 User Manual	
Coll	imat	or	1		
•	Ма	nual Collimator		0232 DX-D	
	0	Push-button to switch on the Collimator Lamp		400 User Manual	
	0	Knobs to open or close the internal blades of the Collimator			
•	Aut	tomatic Collimator (CR)			
	0	Collimator Display			
	0	Motorized Blade Control			
	0	Manual Blade Control			
	0	Automatic Mode Indicator (Green)			
	0	System not ready Indicator (Red)			
	0	Manual Mode Indicator (Yellow)			
	0	Manual Mode Setting of SID			
	0	Cassette-Film Size Memory			
	0	Collimator Lamp Control			
	0	Automatic cassette size sensing – ACSS (if table and wall stand are adapted accordingly)			
Che	ck n	ext page →			

Training Checklist – Manual (Touch Button) Console 7.3.1

Topics: Dor Manual (Touch Button) Console → In case the system is installe softconsole, see next list	11111111111
	d with
SoftConsole, see flext list	
Radiographic and General Controls	0232 DX-D 400 User
o Power on/off	Manual
Workstation Selection	
Focal Spot Indicator	
Radiographic Parameter	
Automatic Exposure Control	
Activation of AEC	
o Exposure Detector Fields (Ion Chamber)	
o Film/Screen Combination	
o Density Compensation	
o AEC Reset	
The Anatomical Programmer (APR)	·
APR display	0232 DX-D 400 User
APR selectors	Manual
APR regions, views and projections	
Patient size	
"Store" button (how to manually rewrite and store the exposure technique)	
"Modify" button (only for key users):	
 Restore Default Replace the data stored in the "APR Working File" with the data stored in the "APR Backup File" 	
 Edit Create, modify or delete the APR techniques Change the operator password Transfer APR Techniques to APR Backup File 	

DX-D 400 - Checklist for Customer Training				
Topics:	Done	Reference		
Exposure Controls and Status Indicators	Exposure Controls and Status Indicators			
X-Ray Handswitch Two stages (Preparation and X-Ray exposure) Advise the radiographer to verify the exposure settings for correctness before doing the exposure.		0232 DX-D 400 User Manual		
"Ready" indicator				
"X-Ray" on				
SELF- DIAGNOSIS INDICATORS				
Door open		0232 DX-D 400		
Generator overload		User Manual		
Tube overload				
Rotor error				
Heat				
• Tech				
End				

Training Checklist - Soft Console 7.3.2

	ftwa fore	are Console → in case system is installed with hardc	onsole s	see checklist
•	0	Operational Status Bar with: - Device status (filtering, grid, alignment status, etc.) - DAP display frame - Error status bar Generator Tab Radiographic Parameters Focal Spot Indicator Automatic Exposure Control - AEC chamber selection - Sensitivity If it is necessary to adapt the cut-off dose, train not to change the Sensitivity, but to use the function Density instead. - Density - Patient thickness		0289 DX-D Software Console User Manual Application Note - NX configuration items in DR systems_max mAs and mA (Document ID: 43122034)
X-ı	ay (Generator Control		
•	E	Buttons to switch on and off		0232 DX-D 400 User Manual
•		(-Ray Handswitch Two stages (Preparation and X-Ray exposure)		Coor Manaci
		Advise the radiographer to verify the exposure settings for correctness before doing the exposure.		
En	d			

7.4 **Detector**

7.4.1 Wireless Detector - DX-D 30C/35C

	DX-D 400 - Checklist for Customer Training Wireless Detector DX-D 30C/35C					
Site	Site: Location: Date:					
Goa Trai	II: ning participant knows how to work with the wireless detect	tor.				
Тор	Topics: Done Reference					
Star	ting the DX-D 30C/35C		0197 EN UM DX-D 30C			
•	Attachment of a fully charged battery pack		DX-D 30C			
•	Switch on the detector					
Stop	oping the DX-D 30C/35C					
•	Turn off detector					
•	Remove battery pack					
Detector status information						
Charging the battery and information about charge status						
	rmation about grid specification and how to attach the grid ional)					
•	Labels on the DR detector		Not			
•	Positioning the DR Detector		applicable			
•	Detector orientation					
•	Charging the detector					
•	Use of protective plastic bag					
•	Cleaning and Disinfection					
End			•			

Wireless Detector - DX-D 40C/G and DX-D 45C/G 7.4.2

	DX-D 400 - Checklist for Customer Training Wireless Detector DX-D 40C/G / DX-D 45C/G				
Site	: Location:	Date:			
Goa Traii	I: ning participant knows how to work with the wireless detect	or.			
Тор	Topics: Done Reference				
Star	ting the DX-D 40C/G or DX-D 45C/G		0290 DX-D 40C.		
•	Attachment of a fully charged battery pack		DX-D 40G		
•	Switch on the detector		User Manual		
Stop	ping the DX-D 40C/G or DX-D 45C/G		0292 DX-D 45C,		
•	Turn off detector		DX-D 45G User Manual		
•	Remove battery pack				
Detector status information					
Cha	Charging the battery and information about charge status				
	mation about grid specification and how to attach the grid onal)				
•	Labels on the DR detector		Not		
•	Positioning the DR Detector		applicable		
•	Detector orientation				
•	Charging the detector				
•	Use of protective plastic bag				
•	Cleaning and Disinfection				
End					

7.4.3 Wireless Detector DR 14C/G and DR10C

DX-D 400 - Checklist for Customer Training Wireless Detector DR 14C/G and DR10C			
Site: Location:	Date:		
Goal:			
Training participant knows how to work with the wireless detect	or.		
Topics:	Done	Reference	
Starting the DR 14C/G and DR10C		0350	
Attachment of a fully charged battery pack		DR 14s Detector	
Switch on the detector		User Manual	
Stopping the DR 14C/G and DR10C			
Turn off detector			
Remove battery pack		0351 DR 10s	
Detector status information		Detector User Manual	
DR detector sharing			
Linking of detector			
Additional necessary information for the customer:			
In case of DR detector sharing where the systems work in different sync mode (Exposync / Auto Trigger Mode): It is required to switch off / on the DR detector after linking.			
Before linking a DR detector on another system, ensure that the last image has been transferred. Linking fails with error messages on NX (e.g. TRIMESS119) if there is still an image on the DR detector. In this case return to the WiFi range of the previous DR detector and wait until the image is transferred.			
Automatic exposure detection			
Charging the battery and information about charge status			
DR 10s / DR 14s specific characteristics to be known by customer:			
When moving the DR 10s / DR 14s DR detector out of the WiFi range and moving back in the WiFi range, it may take up to one minute until the panel is reconnected and the thumbnail is ready for exposure.			
Check next page →			

 For calibration the detector needs to have stable operating temperature. Before each exposure the detector temperature is compared with the temperature during calibration. In case of a deviation of ± 6°C a corresponding warning message is displayed. 	0351 DR 10s Detector User Manual
Stand by (sleep) mode and switch off automatically	
Modes are activated after specific time (configurable)	
Sleep mode can be deactivated with a short delay: After selecting an exposure on the NX workstation.	
Automatically switch off: DR Detector must be started again using the power button.	
Additional necessary information for the customer: The timers do not work in following cases:	
When in NX a thumbnail for DR 10s / DR14s is selected.	
When NX is not running.	
When the DR detector is out of Wi-Fi range.	
During DR detector calibration	
In case of a DR 10s / DR 14s related error	
Information about grid specification and how to attach the grid (optional)	
Labels on the DR detector	Not
Positioning the DR Detector	applicable
Detector orientation	
Charging the detector	
Use of protective plastic bag	
Cleaning and Disinfection	
End	

Checklist for Customer Training "DR Detector Calibration" 7.4.4

Checklist for Customer Training "DR Detector Calibration" DR Detector Calibration				
Site:	Location:	Date:		
Goal: Training	participant is able to calibrate the DR Detector.			
Topics:		Done	Reference	
Calibrati	ng the detector			
Calibratio	on Frequency		0134 DX-D DR Detector	
Calibratio	on Workflow		Calibration Key User Manual	
• Red	quired material			
• Pre	econditions (e.g. warm-up of detector,)			
• Pro	ocedure			
End				

Checklist for Customer Training "NX Workstation" 7.5

Checklist for Customer Training "NX Workstation"				
Site: Location	on:		Date:	
Goal: Training Participant is able to opera	ate the NX work	kstation.		
Topics:			Done	Reference
NX General				
Explanation of the NX User Manual				NX User Manual
Intended use and limitations of the s	system			
Basic Workflow				
Opening a patient from the RIS				NX Key User Manual
Manually entering patient data				
Composing the examination				
Selecting and Performing X-Ray Ex	posures			
DR workflow				
CR workflow				
Performing quality control				1
Editing possibilities E.g. annotations				
"Close and Send" Function			П	
Close and Gend 1 unction				
Worklist Window			ı	
Different panes:				NX User Manual
Search pane				
Worklist pane				
 Closed Exams pane 				
 Manual Worklist pane 				
Image Overview pane				
Action Buttons				
Check next page →				•

Checklist for Customer Training "NX Workstation"		
Topics:		Reference
Browsing through the lists		
Examination Window		
Different panes:		NX User Manual
Patient pane		iviariuar
Image Detail pane		
Image Overview pane		
Action buttons		
Using Examination:		NX User Manual
Image indicators as modality status, status indication of DR detector and problem status		Manual
Tools for performing basic QC (R,L marker, Flip, Rotate)		
Edit patient data		
Edit image data		
Finalizing the exam after the images have been received		
Editing Window		
Advanced editing possibilities:		NX User
Adding annotations to an image and using measurement tools		Manual NX Service
Rotating or flipping an image, show or hide square marker		Manual, chapter 5, "Software
Zooming in or out on an image, applying shutters		Options – Musica²
Working with collimation/cropping, collimation/cropping rules for DR and CR		Image Processing Software"
Cropping (DR) and Black border (CR)		
Use of window/level		
Adjusting the Musica² image processing parameters and explaining the impact of taste adjustments on images		
Check next page →		

Checklist for Customer Training "NX Workstation"		
Topics:		Reference
Editing Window (continued)		
Printing images (including Multi Patient Print)		NX User Manual
Action buttons: Reject, save as new, print sheet, send image, close and send all		
Main Menu		
Different panes:		NX User Manual
Functionality overview pane		iviariuar
Workspace area		
Image Overview Pane		
Action buttons (individually configurable)		
Monitoring and Management		NX User
Queue Management Manua Manua		Iviariuai
Delete Examination		
Lock Examinations		
Quality assurance tasks		NX User
Dose Consistency Reporting		Manual
Repeat Reject Program		
Read and initialize cassettes (optional)		
"Import" and "Export" Functions		
Check next page →		

Checklist for Customer Training "NX Workstation"		
Topics:	Done	Reference
Training for optional Application Packs		
"NX 3.0 Features License":		NX User
Full Screen Mode		Manual
Configurable UI button to call up 3rd party application		
Forced entry of operator ID		
Mandatory Patient fields		
Multiple DICOM export destinations		
Export of technical or clinical images to hard disk		
Print entire study in one click		
Image Quality		
Highlight that the basic rules to create an X-ray exposure are still valid and important for the image quality (e.g. collimation of X-ray beam).		Not applicable
Explain the usage of lead markers.		
Explain the dependencies between X-ray dose, signal-to-noise ratio (SNR) and density.		
Point out that modifications in the image processing may cause image details to get lost (exaggerated MUSICA settings).		
Explain burning/clipping artifacts (exaggerated W/L settings).		
Highlight that the scaling factor of images is independent of real body part size.		Not applicable
Explain the meaning and usage of the "Exposure Index" (EI) as dose indicator, range of over/under exposure.		
Explain the limitation of over/under exposure (dynamic range of a DR detector versus CR): e.g. clipping provoked by a too high exposure		
End		

DX-D 400 Glossary 8

Term	Description
AEC	Automatic Exposure Control
APR	Anatomical Programmer
ACSS	Automatic Cassette Size Sensing
Black border	Function of MUSICA image processing software: The black border masks non-relevant image areas, which are outside the manually collimated image area.
Collimation	Collimation at NX, performed by defining a new region of interest and consecutive re-processing of the image.
	2. Collimation at X-Ray device, made by a diaphragm, an X-ray tube head attachment consisting basically of a sheet of lead with a hole that determines the size and shape of the primary beam.
Cropping	For DR: Region of Interest/Collimated Area is detected with same algorithm as for Black Border CR images. Non-relevant image areas are not masked, but cropped.
DAP	Dose Area Product (DAP), is a multiplication of the dose and the area exposed, often expressed in Gy x cm ² .
meter	Measurement tool to measure the DAP
DI	Deviation Index
DX-D 400 <u>M</u>	M stands for "manual version" of the DX-D 400 system
EI	Exposure Index
Grid	Also called anti-scatter-grid. A device constructed of alternating strips of lead and a radio-transparent medium (such as aluminum or plastic composites) which are oriented in such a way that most of the primary radiation will pass through the grid between the strips while most of the scattered radiation will intersect the lead strips and be absorbed.
Grid ratio	The ratio of the height of the lead grid strips to the distance between the strips.
TEI	Target Exposure Index
μGy	Microgray; an SI unit of absorbed dose equal to 10-6 grays.
	SI = International System of Units
XRDI	X-ray Device Interface

9 AGFA – DX-D 400 - Training and Instruction Protocol

After filling out this form, hand over a copy of the form to the Agfa project manager.

Site Information

[To be adapted locally:

e.g. name of site, name of contact person, site address

Remove this blue text bloc from final version]

DX-D 400 System Description

X-Ray System	CR System
[To be adapted locally: e.g. floor mounted X-ray Unit, Generator,	[To be adapted locally: e.g. digitizer
Remove this blue text bloc from final version]	Remove this blue text bloc from final version]

Objective of Training

DX-D 400 can be operated successfully.

Key Aspects

[To be adapted locally, e.g.:

DX-D 400 workflow

Handling instruction of the X-Ray System

Remove this blue text bloc from final version]

Notes

[To be adapted locally:

Remove this blue text bloc from final version]

AGFA – DX-D 400 Training and Instruction Protocol

Participants List

Date	Name, Title	Signature

Instruction carried out by:		
Date: _		
Signature: _		