

DX-D 100	Type 5410	SB Nr. 46
DX-D 100 Wireless	Type 5411	SB Nr. 43
DX-D 300	Type 8207	SB Nr. 46
DX-D 400	Type 5420	SB Nr. 16
DX-D 600	Type 5430	SB Nr. 36
DX-D 800		SB Nr. 12

DD+DIS048.14E

Service Information Bulletin

This Bulletin is for information only.

Release of Sedecal Techservice 6.15

Task		
Timing		Category
Next service	<input type="radio"/>	Apply at all sites
	<input type="radio"/>	Apply at affected sites as listed below
	<input checked="" type="radio"/>	Optional to improve functionality of product

Task Tracking

After completion of your task the following entry in your Service Report is required:

DD+DIS048.14E

*

* Insert the document number into the field "Comment" (SMS form).

Purpose of this document:

- This document announces the release of version 6.15 of the Sedecal Techservice software.
- It contains information about improved font visibility that has been introduced with Techservice 6.15.

Use, dissemination, distribution or reproduction of this bulletin by unauthorized personnel is not permitted and may be unlawful.
The controlled version of this document is available from the Agfa HealthCare Library. Any printed copy of this document is uncontrolled.

Copyright © 2014 Agfa HealthCare N.V. | Document Node ID: 44987374 | service_information_bulletin_template_v06

Release date: 04-2014 DX-D 100 / Wireless, DX-D 300, DX-D 400,

Agfa HealthCare Company Confidential

DX-D 600, DX-D 800

Page 1 of 4

1 Introduction/purpose

The Sedecal TechService software is used to configure and calibrate the X-ray generator settings of the DX-D 100 / 100 Wireless, DX-D300, DX-D 400, DX-D 600 and DX-D 800 systems..

For configuring the X-ray Generator, install TechService software version 6.15.

For details of the installation steps, refer to section 3 in this document.

1.1 Improved control button visibility on the Techservice GUI

Symptom Currently, the lettering on the controls in the Techservice GUI is almost invisible after disabling a control's function.

Cause The coloring of the GUI items.

Solution In Techservice 6.15, the color of the lettering has been changed so that it remains visible when a control's function is disabled.

2 Prerequisites



SOFTWARE:

Download the **TechServiceInstall6.15.exe** software (file size 2399 KB) from the *AGFA HealthCare Library* to a portable storage medium:

Direct Radiography > DX-D Components > Generator > Software



TOOLS:

- Portable storage medium (checked to be virus-free, e.g. USB Memory Stick)



REFERENCED DOCUMENTS:

Download the documents from the *AGFA HealthCare Library*:

- DX-D 100 / DX-D 100 Wireless - Mobile X-Ray Unit - Service Manual
Direct Radiography > DX-D 100 > Service Manual
- X-ray Generator – Service Manual – DX-D 300
- X-ray Generator – Service Manual – DX-D 400
- X-ray Generator – Service Manual – DX-D 600
- X-ray Generator – Service Manual – DX-D 800
Direct Radiography > DX-D Components > Generator > Service Manual

3 Instructions



REQUIRED TIME:

Approximately 30 minutes



NOTE:

The Techservice software can be installed by running the executable. In case a version of the Techservice software is already installed, please remove this first before installing the new version.

- (1) Stop the NX application before installing Techservice version 6.15.
Go to: **Start > Agfa > NX > Stop NX.**
- (2) **Click Software TechServ 6.15.exe** and follow the instructions on the screen.
It is possible that a popup will appear after installation. If so, select the option **This program installed correctly.**



4 Verification

- (1) Restart the NX application:
Start > Agfa > NX > Restart NX Completely
- (2) Go to **Start > Settings > Control Panel > Add** or **Remove Programs** and verify that the installed Techservice software version is 6.15

5 Keywords

Techserv, Techservice 6.15, TechServiceInstall6.15.exe

6 Version history

Version	Change	Date
1.0	Initial version	04-2014