

<b>DX-D 100</b>	Type 5410	SB No. 31
<b>DX-D 100 Wireless</b>	Type 5411	SB No. 28
<b>DX-D 300</b>	Type 8207	SB No. 37
<b>DX-D 400</b>	Type 5420	SB No. 08
<b>DX-D 500</b>	Type 8206/050	SB No. 29
<b>DX-D 600</b>	Type 5430	SB No. 22
DD+DIS102.13E		

## Service Bulletin

### Announcement of spare part silicon paste for maintenance

Task			
Timing		Category	Scope
Next service as agreed with customer	<input checked="" type="radio"/>	Apply at all sites	PowerHelp complaint / HQ issue: HQ_1211280002
	<input type="radio"/>	Apply at affected sites as listed below	
	<input type="radio"/>	Optional to improve functionality of product	

### Task Tracking

After completion of your task the following entry in your Service Report is required:

DD+DIS102.13E

\*

\* Insert the document number into the field "Comment" (SMS form).

Purpose of this document:

This document announces the spare part silicon paste for maintenance tasks.

The signatures on the approval page indicate the solutions described in this Service Bulletin have been reviewed and are NOT reportable because no actions are taken to reduce a "Risk to Health" according to our risk assessment process.

## 1 Introduction/purpose

<i>Symptom</i>	It is not possible to power the DX-D System on after maintenance actions at the X-Ray tube.
<i>Cause</i>	During maintenance of the X-Ray tube the HV cables will be disconnected. When reconnecting them it is possible, that air bubbles occur at the high tension cables and the connection is not fully established. As a result the DX-D System can't be powered on.
<i>Solution</i>	Use the silicon paste to lubricate the connectors to avoid air bubbles. That guarantees full contact of the connectors and avoids damage to the connectors. Presence of air bubbles can cause electrical arcing which can cause damage to the HV cables.

## 2 Prerequisites



SPARE PARTS:  
CM+6074010 Silicon Paste for X-Ray tubes



TOOLS:

- Eye protection
- Protective clothes

### 3 Instructions



**REQUIRED TIME:**  
5 minutes



**NOTE:**

Always use the silicone paste after a maintenance action at the X-Ray tube.

- (1) Switch off the DX-D System.
- (2) Disconnect mains cable.
- (3) Perform maintenance as described in the corresponding manual.
- (4) When reconnecting the HV cables:  
Apply the silicone paste on the connectors of the HV cables at the X-Ray tube before reconnecting them.
- (5) Reconnect mains cable.
- (6) Switch on the DX-D System.  
If the system does not start, refer to other system relevant service bulletins and documentation.

### 4 Verification

- (1) Power on the system.

### 5 Keywords

Silicon paste, HV cable, high tension, X-Ray tube

### 6 Version history

Version	Change	Date
1.0	Initial version	10-2013



## Details as of PDF Creation Date

### Document Metadata

<b>Title:</b>	DX-D 100 - Service Bulletin No. 31 - Announcement of spare part silicon paste for maintenance
<b>Livelink ID:</b>	42204272
<b>Version#:</b>	9
<b>Version Date:</b>	2013-10-08 01:55 PM CET
<b>Status:</b>	Approved on 2014-03-12 01:19 PM CET
<b>Owner:</b>	Carsten Bockel (awyuq)
<b>Created By:</b>	Carsten Bockel (awyuq)
<b>Created Date:</b>	2013-09-06 04:08 PM CET
<b>PDF Creation Date:</b>	2014-03-12 01:19 PM CET

**This document was approved by:**

### Signatures:

1. Josef Wagner (agwj) on 2013-10-08 04:17 PM CET
2. Paul Merckx (amdag) on 2013-10-09 08:51 AM CET
3. Bart Biesemans (amajm) on 2013-10-10 09:58 AM CET
4. Lieven Lauwers (awibr) on 2013-10-15 11:07 AM CET
5. Christian HÃ¼hne (aghch) on 2014-03-12 10:58 AM CET

### Detailed Approver History:

- **Approval Workflow started on 2013-10-08 02:53 PM CET**
  - Approval task originally assigned to and completed by Josef Wagner (agwj) on 2013-10-08 04:17 PM CET
  - Approval task originally assigned to and completed by Christian HÃ¼hne (aghch) on 2014-03-12 10:58 AM CET
  - Approval task originally assigned to and completed by Lieven Lauwers (awibr) on 2013-10-15 11:07 AM CET
  - Approval task originally assigned to and completed by Paul Merckx (amdag) on 2013-10-09 08:51 AM CET
  - Approval task originally assigned to and completed by Bart Biesemans (amajm) on 2013-10-10 09:58 AM CET

### Version & Status History

--	--	--

Version#	Date Created	Status
9	2013-10-08 01:55 PM CET	Approved - 2014-03-12
8	2013-10-08 01:54 PM CET	
7	2013-10-07 08:39 AM CET	Reviewed - 2013-10-08
6	2013-10-01 09:24 AM CET	
5	2013-10-01 09:14 AM CET	
4	2013-09-13 03:37 PM CET	Reviewed - 2013-09-26
3	2013-09-13 03:18 PM CET	
2	2013-09-06 04:09 PM CET	
1	2013-09-06 04:08 PM CET	