

## Service Information Bulletin

This Bulletin is for information only.

### Driver for first version of WiFi cards available to solve sporadic blue screen after clean installation

#### Task

Timing		Category	
Next service	<input type="radio"/>	Apply at all sites	
	<input checked="" type="radio"/>	Apply at affected sites as listed below	
	<input type="radio"/>	Optional to improve functionality of product	

#### Task Tracking

After completion of your task the following entry in your Service Report is required:

DD+DIS145.13E

\*

\* Insert the document number into the field "Comment" (SMS form).

Purpose of this document:

- This document informs about a sporadic "blue screen" at devices with "Via Technologies Wifi Card".
- Solution: Install the WiFi card driver which is available on the Agfa HealthCare Library.

Affected serial number(s) / batch:

All DX-D 100 / DX-D 100 Wireless manufactured until 03-2012 may be affected, where a clean installation with Windows 7 has been performed.

## 1 Introduction

This document informs about the solution for the following software issue:

**Symptom** Sporadically after startup or re-login after a screen lock, a “blue screen” is shown. The mouse pointer still works but no other programs are accessible. The system is unusable. A reboot via main switch is required.

Both of the following prerequisites must be fulfilled for the symptom to appear:

- DX-D 100 or DX-D 100 Wireless has been installed from scratch in the field via Windows 7 clean installation DVD (DXD100 Base OS image1.1\_Win7.iso).
- DX-D 100 or DX-D 100 Wireless has the “Via technologies WiFi” card installed. See also NOTE below.



**NOTE:**

Two types of WiFi Lan cards for communication with the hospital network have been installed in production:

- Until 03-2012: Via Technologies Wifi Card
- As of 03-2012: Sparklan Wifi Card (= 802.11n USB Wireless LAN Adapter)

These cards require different drivers.




Figure 1

**Cause** The Windows 7 operating system provided for DX-D 100 / DX-D 100 Wireless on the Fileshare FTP server has USB WLAN driver VT6656\_Win7\_V1.1.0.2\_32bit pre-installed. This driver may cause sporadic “blue screens”.

**Solution** Install the previous version of this driver: VT6656 V1.0.0.2.  
For details see next sections.

## 2 Prerequisites



### SOFTWARE:

Download the DX-D 100 USB wireless Lan driver VT 1.0.0.2 from the Agfa HealthCare Library:

Direct Radiography > DX-D 100 Wireless > DX-D 100 Wireless Generic or

Direct Radiography > DX-D 100 > DX-D 100 Generic



### TOOLS:

- USB flash drive
- USB mouse

## 3 Instructions



### REQUIRED TIME:

Approx. 30 minutes:

- Approx. 10 minutes installation of the driver
  - Approx. 20 minutes verification
- (1) Check WiFi communication, e.g. by checking the Network connections in the explorer. This allows comparing the WiFi communication before and after the driver downgrade.
  - (2) Go to **Control Panel > Device Manager > Network Adapters**.
  - (3) Check the type of Wireless LAN card: Only proceed if it is a Via Technologies Adapter. See also Figure 1 on previous page.
  - (4) Right-click on **Via Technologies USB Wireless LAN adapter**.
  - (5) Click **Properties**.
  - (6) Go to the Driver tab: The driver version must be 1.0.0.2.




Figure 2

- (7) If a different driver version is displayed, right-click on **Via technologies USB Wireless LAN adapter** and select **Uninstall**.
- (8) Unzip the DX-D 100 USB WLAN driver VT 1.0.0.2 and copy file "VT6656\_Win7\_V1.0.0.2\_32bit.exe" via USB flash drive to the NX.
- (9) Execute the driver.
- (10) In case the Program Compatibility Assistant appears, select **This program installed correctly**.
- (11) Again check the driver version as described above: It must be 1.0.0.2.

## 4 Verification

- (1) Check WiFi communication, e.g. by checking the network connections in the explorer.
- (2) Stop NX.
- (3) Reboot the PC.
- (4) Check that no "blue screen" appears after reboot.
- (5) Repeat steps (2) to (4) a few times.

## 5 Keywords

Bluescreen, clean install

## 6 Version history

Version	Change	Date
1.0	Initial version	09-2013