

## Service Bulletin

### Battery Overheating Issue

Task			
Timing		Category	Scope
Next service as agreed with customer	<input checked="" type="radio"/>	Apply at all sites	PowerHelp complaint / HQ issue: HQ_1204260002
	<input type="radio"/>	Apply at affected sites as listed below	
	<input type="radio"/>	Optional to improve functionality of product	

#### Task Tracking

After completion of your task the following entry in your Service Report is required:

DD+DIS123.12E

\*

\* Insert the document number into the field "Comment" (SMS form).

Purpose of this document:

- To inform about the addition of a thermostat in the DX-D 100 / DX-D 100 Wireless that prevents overheating of a battery and the adjacent batteries during the charging process.
- To inform about availability of the corresponding new schematics.

The signatures on the approval page indicate the solutions described in this Service Bulletin have been reviewed and are NOT reportable because no actions are taken to reduce a "Risk to Health" according to our risk assessment process.

Use, dissemination, distribution or reproduction of this bulletin by unauthorized personnel is not permitted and may be unlawful.  
The controlled version of this document resides on the HealthCare Library. Any printed copy of this document is uncontrolled.

Copyright © 2012 Agfa HealthCare N.V. | Document Node ID: 36874631 | service\_bulletin\_template\_v06

## 1 Introduction/purpose

*Symptom* Sulfation (forming of lead sulfate powder) and overheating of batteries next to the battery which is being charged, caused by heat which is transmitted to the other batteries.

*Cause* When the battery of the DX-D 100 / DX-D 100 Wireless is charging for a very long time there is a possibility that it will overheat during the charging process.

*Solution* A thermostat will be added inside the DX-D100 / DX-D 100 Wireless.  
This thermostat will interrupt the charge of the batteries if the temperature inside the DX-D 100 reaches 50°C (122° F).  
This prevents heat damage to the charged battery and the batteries next to it.

### NOTE:



The newer DX-D 100 / DX-D 100 Wireless solutions are delivered with a thermostat.  
Older DX-D 100 / DX-D 100 Wireless solutions which were delivered before a certain point in time may not have been fitted with a thermostat.



### IMPORTANT:

So, check whether a thermostat has been installed in the solution or not.  
If no thermostat has been installed, follow the procedure described below.

## 2 Prerequisites



### SPARE PARTS:

The Thermostat Service Kit can be ordered via following spare part number\*:  
SC+A8224-01

- \* The last digit in the spare part number indicates the spare part revision at release of this document.  
When ordering, the actual revision of the spare part is delivered.



### REFERENCED DOCUMENTS:

- DX-D 100 Wireless schematics:  
54302150, revision I (refer to document DX-D 100 Wireless - X-Ray Mobile –  
Schematics in the Agfa HealthCare Library):  
**Direct Radiography → DR Equipment → DX-D 100 Wireless**
- DX-D 100 schematics:  
54302200, revision E (refer to document DX-D 100 - X-Ray Mobile –  
Schematics in the Agfa HealthCare Library):  
**Direct Radiography → DR Equipment → DX-D 100**

### 3 Instructions

**REQUIRED TIME:**

Approximately 45 minutes for mounting

#### 3.1 Mechanical Preparations at the Console

- (1) Disconnect the device from the mains.
- (2) Remove console frame:
  - Open the pencil tray and remove the two screws that fasten the console frame to the unit.
  - Push the console frame towards the column. Do not lift up the console frame until its upper side is fully released.



Figure 1

- (3) Remove the front covers:
  - Remove the six lateral screws placed on the front cover (see dotted-line marking) at both sides underneath the console frame. See figures 2 and 3.



Figure 2



Figure 3

- (4) Pull the curved-top cover upwards to disengage the three snaps underneath.



Figure 4

- (5) Remove the two fixing screws that fix the front cover to the mobile frame and remove the three screws to dismount the Bucky connector cover.




Figure 5

- (6) Before dismounting the front cover, remove two screws to dismount the catch of the parking detent device; then disconnect the sensor signal cable connector (J13).



Figure 6

- (7) Lift the front cover from the unit as shown in figure 7, passing the power line plug through its window. For units with infrared remote control (optional), disconnect the sensor signal cable from the PCB and pass it through the cover hole.



Figure 7

- (8) Loosen the screws of the LF1 line filter behind which the thermostat will be fixed so that the thermostat support can be slid behind the line filter.
- (9) Dismount thermostat covers.

The following 2 figures illustrate how the thermostat kit is installed behind the line filter:

- (10) Mount the thermostat on the support which is fixed at the mobile location shown in figure 10 (dotted red line).




Figure 8

- (11) Figure 11 shows the position of the bracket (with the installed switch) mounted behind the line filter in greater detail. The area marked with a dotted frame in figure 12 corresponds exactly with the detail in figure 11.




Figure 9




Figure 10

## 3.2 Electrical Installation

The following two figures show the complete wiring diagram (fig. 13) and the relevant detail which shows the new wiring layout with the built-in thermostat (figure 14).




Figure 11




Figure 12





## 4 Verification

Perform a complete battery charging cycle and check if the issue occurs again.

The thermostat should interrupt the charging process if the temperature inside the DX-D 100 reaches 50°C (122° F).

## 5 Keywords

Battery, overcharge, thermostat, charging, interrupt.

## 6 Version history

Version	Change	Date
1.0	Initial Version	07-2012



## Details as of PDF Creation Date

### Document Metadata

<b>Title:</b>	DX-D 100 Wireless - Service Bulletin No. 11 - Battery Overheating Issue
<b>Livelink ID:</b>	36874631
<b>Version#:</b>	14
<b>Version Date:</b>	2012/07/13 02:04 PM CET
<b>Status:</b>	Approved on 2012/07/16 02:38 PM CET
<b>Owner:</b>	Beate Richter (axnwp)
<b>Created By:</b>	Beate Richter (axnwp)
<b>Created Date:</b>	2012/05/02 03:37 PM CET
<b>PDF Creation Date:</b>	2012/07/16 02:38 PM CET

**This document was approved by:**

### Signatures:

1. Josef Wagner (agwj) on 2012/07/13 02:48 PM CET
2. Bart Biesemans (amajm) on 2012/07/13 02:44 PM CET
3. Pieter Perdieu (amaft) on 2012/07/13 04:14 PM CET
4. Lieven Lauwers (awibr) on 2012/07/16 02:35 PM CET

### Detailed Approver History:

- **Approval Workflow started on 2012/07/13 02:14 PM CET**
  - Approval task originally assigned to Paul Merckx (amdag)
    - Approval task completed by Pieter Perdieu (amaft) on 2012/07/13 04:14 PM CET (workflow proxy)
  - Approval task originally assigned to and completed by Lieven Lauwers (awibr) on 2012/07/16 02:35 PM CET
  - Approval task originally assigned to and completed by Josef Wagner (agwj) on 2012/07/13 02:48 PM CET
  - Approval task originally assigned to and completed by Bart Biesemans (amajm) on 2012/07/13 02:44 PM CET

### Version & Status History

Version#	Date Created	Status

14	2012/07/13 02:04 PM CET	Approved - 2012/07/16
13	2012/07/13 02:00 PM CET	
12	2012/07/11 04:53 PM CET	Approved - 2012/07/13
11	2012/07/11 04:48 PM CET	
10	2012/06/18 03:35 PM CET	Approval Cancelled - 2012/06/29
9	2012/06/18 03:34 PM CET	
8	2012/05/31 03:23 PM CET	Approved - 2012/06/11
7	2012/05/31 02:50 PM CET	
6	2012/05/23 02:16 PM CET	
5	2012/05/23 01:59 PM CET	