

## Service Bulletin

### Default Exam Tree for DX-D 100 available

Task			
Timing		Category	Scope
Next service as agreed with customer	<input checked="" type="radio"/>	Apply at all sites	PowerHelp complaint / HQ issue: HQ_1010080002
	<input type="radio"/>	Apply at affected sites as listed below	
	<input type="radio"/>	Optional to improve functionality of product	

### Task Tracking

After completion of your task the following entry in your Service Report is required:

DD+DIS343.10E

\*

\* Insert the document number into the field "Comment" (SMS form).

Purpose of this document:

This Service Bulletin communicates the availability of a default exam tree for every DR Solution.

The signatures on the approval page indicate the solutions described in this Service Bulletin have been reviewed and are NOT reportable because no actions are taken to reduce a "Risk to Health" according to our risk assessment process.

Use, dissemination, distribution or reproduction of this bulletin by unauthorized personnel is not permitted and may be unlawful.  
The controlled version of this document resides on MedNet. Any printed copy of this document is uncontrolled.

Copyright © 2010 Agfa HealthCare N.V. | Document Node ID: 31802106 | service\_bulletin\_template\_v06

## 1 Introduction/purpose

- Symptom* When importing a DR exam tree of another DR Solution, the density steps of some exposures get automatically (wrongly) set to density value -3.
- Cause* Mixing the exam trees in the different DR solutions can cause changes in the density steps.
- Solution* A default exam tree for every DR solution is available.



### IMPORTANT:

Make sure only to use the DR solution specific exam tree.  
It is not allowed to use an exam tree designed for another DR Solution.

## 2 Prerequisites



### TOOLS:

Portable storage medium (checked to be virus-free, e.g. USB Memory Stick)



### SOFTWARE:

Download the default exam tree (File size approximately 2 MB) from the MedNet GSO Library Path to a portable storage medium ([Intranet Link](#) / [Extranet Link](#)).

Direct Radiography => DR Equipment => DX-D 100 => Software => DX-D 100 - Default Exam tree



### REFERENCED DOCUMENTS:

Application System Manual - Chapter 6 - Configuration & Customization

### 3 Instructions

**REQUIRED TIME:**

Download time exam tree : 5 min

- (1) Open the NX Configuration Tool.
- (2) Go to: **Exam Tree Configuration**
- (3) Select: **<Add Exams>** from Import...
- (4) Browse to the exam tree stored on your portable storage medium
- (5) Open the default Exam Tree from that location.
- (6) Click **<OK>**.
- (7) Configure and Customize the Exam Tree as described in Application System Manual - Chapter 6 - Configuration & Customization.
- (8) Select *File → Activate configuration* in the configuration tool menu.

### 4 Verification

**With activation of the configuration the configuration is automatically verified.**

- (1) Solve all configuration errors reported by the Configuration Tool (if there are some).
- (2) Retrigger the activation of the configuration.

### 5 Keywords

Exam tree, density value

### 6 Version history

Version	Change	Date
1.0	Initial Version	12-2010



## Details as of PDF Creation Date

### Document Metadata

<b>Title:</b>	SIB02_DD+DIS343.10E_DX-D100 Default Exam Tree
<b>Livelink ID:</b>	31802106
<b>Version#:</b>	9
<b>Version Date:</b>	2010/12/09 02:34 PM CET
<b>Status:</b>	Approved on 2010/12/10 01:57 PM CET
<b>Owner:</b>	Michael Roeser (exfrm)
<b>Created By:</b>	Michael Roeser (exfrm)
<b>Created Date:</b>	2010/11/19 09:13 AM CET
<b>PDF Creation Date:</b>	2010/12/10 01:58 PM CET

**This document was approved by:**

### Signatures:

1. Josef Wagner (agwj) on 2010/12/09 04:53 PM CET
2. Paul Merckx (amdag) on 2010/12/10 10:25 AM CET
3. Bart Biesemans (amajm) on 2010/12/09 04:16 PM CET
4. Frederik Cornelis (amdab) on 2010/12/10 07:28 AM CET

### Detailed Approver History:

- **Approval Workflow started on 2010/12/09 02:36 PM CET**
  - Approval task originally assigned to and completed by Bart Biesemans (amajm) on 2010/12/09 04:16 PM CET
  - Approval task originally assigned to and completed by Josef Wagner (agwj) on 2010/12/09 04:53 PM CET
  - Approval task originally assigned to and completed by Frederik Cornelis (amdab) on 2010/12/10 07:28 AM CET
  - Approval task originally assigned to and completed by Paul Merckx (amdag) on 2010/12/10 10:25 AM CET

### Version & Status History

Version#	Date Created	Status
9	2010/12/09 02:34 PM CET	Approved - 2010/12/10
8	2010/12/09 02:34 PM CET	

7	2010/12/03 09:22 AM CET	
6	2010/12/03 09:11 AM CET	
5	2010/11/23 01:12 PM CET	Reviewed - 2010/12/03
4	2010/11/19 09:18 AM CET	
3	2010/11/19 09:17 AM CET	
2	2010/11/19 09:14 AM CET	
1	2010/11/19 09:13 AM CET	