

Service Information Bulletin

This Bulletin is for information only.

Intermittent failure: Exposure parameters do not appear within the GUI of the examination window.

Task		
Timing		Category
Next service	<input type="radio"/>	Apply at all sites
	<input checked="" type="radio"/>	Apply at affected sites as listed below
	<input type="radio"/>	Optional to improve functionality of product

Purpose of this document:

It informs about an effect that appears in case the "Edit Image" functionality is used by the operator prior to receiving the exposure parameters.

Affected serial number(s) / batch: All DX-D100 in combination with NX 2.0.8400 / 3.0.8400

Use, dissemination, distribution or reproduction of this bulletin by unauthorized personnel is not permitted and may be unlawful.
The controlled version of this document resides on MedNet. Any printed copy of this document is uncontrolled.

Copyright © 2011 Agfa HealthCare N.V. | Document Node ID: 33032778 | service_information_bulletin_template_v05

1 Introduction/purpose

Symptom Exposure parameters do not appear within the GUI of the examination window.

Cause If a patient is selected from RIS or manually entered the required routine is used for the exam thumbnail to appear, i.e Chest PA.

If the Edit Image functionality ("Edit" function within the 'Examination' window) is used to change the original exam setting prior to receiving the exposure parameters, the exposure parameters do not appear within the GUI.

On edit exposure the exposure parameters are not resend.

XRDI still replies with the original SubExposureUid which is not found anymore in the active presentation state. Consequently the received parameters are not stored anymore.

Solution Will be solved in NX 2.0.8600 / 3.0.8600.



IMPORTANT:

Inform the customer:

The "Edit" function within the "Examination" window should not be used prior the exposure parameters are shown within the GUI of the examination window.

2 Prerequisites

Not applicable.

3 Instructions

Not applicable.

4 Verification

Not applicable.

5 Keywords

Exposure parameters, GUI of the examination window

6 Version history

Version	Change	Date
1.0	Initial Version	07-2011