

## Service Information Bulletin

This Bulletin is for information only.

### Sellable upgrades for DX-D 100 / DX-D 100 Wireless available

Task		
Timing		Category
Next service	<input type="radio"/>	Apply at all sites
	<input type="radio"/>	Apply at affected sites as listed below
	<input checked="" type="radio"/>	Optional to improve functionality of product

Purpose of this document:

- This document contains all necessary information about sellable upgrade kits for the installed base DX-D 100 / DX-D 100 Wireless.
- It lists the corresponding upgrade instructions.

## 1 Introduction/purpose

With the recent releases for DX-D 100 / DX-D 100 Wireless quite a number of technical improvements and features have been introduced.

Some upgrade kits have been created to provide access for interested customers to the latest evolutions for the installed base.

These upgrades consist of a combination of "spare part" + working time (if applicable).

The upgrade kits are available both through the Sales Channel (as an ABC Code) and the Service Channel (as a spare part number).

As these upgrade kits are no spare parts as such, the spare part numbers in section 2 will not be included in the spare parts list.

*Supporting material*

The attached leaflet is available to promote and explain the upgrades to the customer:



DX-D 100\_Upgrade  
leaflet.pdf

The leaflet is also available from the Agfa HealthCare Library, Document ID [51330970](#)

## 2 Prerequisites



### PARTS:

The following table contains the order numbers of the upgrades and lists the corresponding installation instructions.

For document details see "Referenced Documents".

Order number	Description	Installation instructions
<b>10+9999913110</b>	DX-D 100 Wireless only: <b>DX-D 100 upgrade to redesigned storage bin</b> The new storage bin is more robust and leaves more space for the detectors.	SB 54 - Announcement of DX-D 100 Version 5.0
<b>10+9999913120</b>	<b>DX-D 100 Upgr. collimator dose optim.</b> Collimator (Ralco R221F) with Double Laser for Optical SID	SB 39 - Upgrade of DX-D 100 / DX-D 100 Wireless to System Version 4.4
<b>10+9999913130</b> (+ Generator SN) *	<b>DX-D 100 Upgr. DAP meter</b> Kit for installing a Vacudap 1580013 internal dose meter	
<b>10+9999913140</b>	<b>DX-D 100 Upgr. IR remote</b> Remote Exposure Control Kit (Infrared wireless remote control)	
<b>10+9999913150</b>	<b>DX-D 100 Upgr. Crystal batteries</b> Upgrade Kit for generator Crystal lead battery	SB 68 - Replacement instructions for generator batteries
<b>10+9999913160</b>	<b>Magnetic holder Sani wipe</b>	SB 54 - Announcement of DX-D 100 Version 5.0
<b>10+9999913170</b>	<b>Magnetic holder hand sanitizer</b>	
<b>10+9999913180</b>	<b>Magnetic holder gloves</b>	

\* Provide the generator serial number together with the order.



### NOTE:

Delivery times may be longer than for regular spare parts.

**REFERENCED DOCUMENTS:**

The documents are available in the Agfa HealthCare Library:

- DX-D 100 - Service Bulletin No. 68 - Replacement instructions for generator batteries, DD+DIS087.15E, Document ID [49961765](#).
- DX-D 100 - Service Bulletin No. 54 - REVISED - Announcement of DX-D 100 Version 5.0, Document ID [45533923](#).
- DX-D 100 - Service Bulletin No. 39 - Upgrade of DX-D 100 / DX-D 100 Wireless to System Version 4.4, DD+DIS006.14E, Document ID [43910047](#).

## 3 Instructions

### 3.1 Inform the Customer

Inform the customer about the availability of the upgrades. You can give the leaflet to the customer to explain the upgrades.

### 3.2 Upgrade installations

In case of required upgrades:

Order the upgrades with the order numbers listed in section 2.

The installation is described in the corresponding installation instructions, listed in section 2.

## 4 Verification

For verification instruction refer to the referenced documents.

## 5 Keywords

Sellable upgrade, optional, service channel

## 6 Version history

Version	Change	Date
1.0	Initial Version	02-2016