

Service Information Bulletin

This Bulletin is for information only.

Announcement of up to date DR Detector related Information stored on FTP Server

Task		
Timing		Category
Next service	<input type="radio"/>	Apply at all sites
	<input type="radio"/>	Apply at affected sites as listed below
	<input checked="" type="radio"/>	Optional to improve functionality of product

Task Tracking

After completion of your task the following entry in your Service Report is required:

DD+DIS279.10E

*

* Insert the document number into the field "Comment" (SMS form).

Purpose of this document:

This document informs about an intermediate approach how to provide up to date information of the DR Detector in a fast way.

Use, dissemination, distribution or reproduction of this bulletin by unauthorized personnel is not permitted and may be unlawful.
The controlled version of this document resides on MedNet. Any printed copy of this document is uncontrolled.

Copyright © 2010 Agfa HealthCare N.V. | Document Node ID: 31214742 | service_information_bulletin_template_v05

1 Introduction

**NOTE:**

The DR Detector Service Manual for the DX-D 100 is included in the appropriate chapter of the DX-D 100 Service Documentation.

In addition to the DX-D 100 Service Documentation R&D representatives provide additional DR Detector information. An intermediate approach has been set up to distribute this information via a FTP Server partition.

For details how to obtain this information via the FTP Server, see section 2.

2 FTP Server Access Prerequisites

**PREREQUISITES:**

- Broadband internet connection
- FTP client like CuteFTP, WSFTP, Filezilla, SmartFTP or any similar application

3 FTP Server Access

- (1) Open your FTP client.
- (2) Connect to the FTP server using following data:

Host	ftp.agfa.be
User	Varian-10x
Password	Detector
Port	21

- (3) Check for latest DR Detector information.
This is indicated via <Date>.
- (4) Select and download file.
- (5) Be aware of the provided DR Detector description.

**IMPORTANT:**

If a connection via FTP client cannot be established try different 'passive transfer mode' settings for your FTP client.

**NOTE:**

For Agfa employees only:

From outside Europe there can be connection problems using the Agfa VPN.

Do **not** download files over the Agfa VPN. Instead connect the FTP server directly from the internet.

**NOTE:**

Check regularly the FTP Server for latest information.

4 Keywords

- DR Detector
- FTP Server

5 Version History

Version	Change	Date
1.0	Initial Version	10-2010