

Radiology Solutions Services

DX-D 100

**DX-D 100 Wireless** 

Type 5410

No. 128

Type 5411 No. 166

Document ID 75129232

# **Service Information Bulletin**

This Bulletin is for information only.

# New spare parts "Hard disk drive" and "Power Supply" for the integrated NX PC available

#### **Task**

Timing			Category
	Next service		Apply at all sites
		•	Apply at affected sites as listed below
		0	Optional to improve functionality of product

## **Task Tracking**

After completion of your task the following entry in your Service Report is required:	75129232	*
* Insert the document number into the field "Comment" (SMS form).		

## Purpose of this document:

This document contains all necessary information about the introduction of new spare parts for DX-D 100 / DX-D 100 Wireless integrated NX PC.

## Affected systems:

All systems where the hard disk drive or power supply of the integrated Windows 7 PC or Windows 10 PC is defective.

Use, dissemination, distribution or reproduction of this bulletin by unauthorized personnel is not permitted and may be unlawful. The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

# 1 Introduction

Reduce the order of DX-D 100 integrated NX PC in case of problems. Supplier's repair analysis of the DX-D 100 NX PC showed that in many cases a replacement of the PCs HDD or the PCs Power Supply could have solved the problem. Since these parts only have fraction of the cost of a complete PC, it was decided to offer the following components as separate spare parts from now on.

New spare part	Compatible with	Picture
Power Supply Computer Mobile SC+53418076	Windows 7 and Windows 10 PC in DX-D 100 / DX-D 100 Wireless	
SSHD 1Tb/8Gb 2,5", SATA III (WIN7) SC+55001313	Windows 7 PC in DX-D 100 and DX-D 100 Wireless. See also NOTE below.	2.3.2.2.2.2.2.2.2.2.2.2.2.2.2.2.2.2.2.2
SSHD 1Tb/8Gb 2,5", SATA III (WIN10) SC+55001403	Windows 10 PC in DX-D 100 Wireless. See also NOTE below.	A STATE OF THE PARTY OF THE PAR



#### NOTE:

- DX-D 100: Was delivered with Windows 7 only (initially Windows XP).
- DX-D 100 Wireless: Was initially delivered with Windows 7.
   Windows 10 is installed for SN ≥ A5411003212 and additionally in A5411003210



#### **REFERENCED DOCUMENTS:**

- DX-D 100 Wireless V11 Service Manual Document ID 66210593
- DX-D 100 V11 Service Manual Document ID <u>66209499</u>
- DX-D 100 / DX-D 100 WL-Mobile X-Ray Unit Service Manual, Document ID 41600560
- DX-D 100 / DX-D 100 WL-Mobile X-Ray Unit Spare Parts List, Document ID 30872155

# 2 Verification

Not Applicable

# 3 Keywords

Workstation, Spare Parts, NX PC

# 4 Version history

Version	Change	Date
1.0	Initial version	12-2020