


Radiology Solutions Services

Document ID: 31372837

DX-D Portable DR Detector – Part Number 7358 **DX-D Portable DR Detector – Part Number 20665**

Type 5410/100



DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.



Manufacturer

Agfa NV

Published by

Agfa-Gevaert HealthCare GmbH
Tegernseer Landstraße 161
D - 81539 München
Germany

Copyright © 2022 Agfa NV

All rights reserved.

Technical modifications reserved.

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert N.V., Belgium, or its affiliates.

All other trademarks mentioned in this document are held by Agfa NV or the respective owners and are used in an editorial fashion with no intention of infringement.

Nothing contained in this legal notice nor in any text in this document shall be construed as granting by implication, estoppel or otherwise, any license or right to use any of the trademarks, service marks, trade names or logos appearing in this document without the express prior written consent of their respective owner.



WARNING:

Improper operation or service activities may cause damage or injuries.

- (1) Read the *Generic Safety Directions* prior to attempting any operation, repair or maintenance task on the equipment. Refer to Document ID 11849633.
- (2) Strictly observe all safety directions within the *Generic Safety Directions* and on the product.



WARNING:

Hazards may be introduced because of component failure or improper operation.

- (1) Replace defective parts with Agfa original spare parts.
- (2) Use only tools and measuring instruments which are suitable for the procedure.
- (3) Only approved Agfa accessories must be used. For a list of compatible accessories contact your local Agfa NV organization or www.agfa.com.

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.



Changes compared to previous revision

The following modifications have been implemented:

- Spare Parts added



DISCLAIMER:

Actual delivered spare parts may differ (in appearance) from the images shown and/or ordered spare parts.

The codes of the ordered spare parts can vary from the code of the delivered spare parts.

Agfa NV warrants that the delivered spare parts have at least the same functionalities.

Agfa NV reserves the right to deliver compatible or alternative spare parts.

Prices may differ from the original order as these spare parts will be invoiced at current prices.



NOTE:

Repaired Spare Parts are indicated with Prefix RF+.

Refurbished Spare Parts are indicated with Prefix CF+.



NOTE:

Recycling of the electronic and electrical waste equipment will ensure safety of the human health and the environment.

For information about electronic and electrical waste equipment disposal, recovery and collection points, please contact your local waste disposal service or the producer / distributor of this equipment.

If your equipment contains removable batteries or accumulators please dispose of these separately according to local regulations.



Contact

Spare Parts ordering


Worldwide: spareparts@agfa.com

Spare Parts returns

Worldwide: returnsspareparts@agfa.com

DOCUMENT CONTROL NOTE:


The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.



01



02



03



04



05




06



07

07, 09

08, 10, 11



12



MOBILE X-RAY UNIT

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

Item No.	Part No.	Description	Replacement Instruction ID
1	RF+29888	DX-D 10G CASSETTE DETECTOR, RECEPTOR DRZ+ 4336R REV. A AND B - (REP.)	-
1	VA+29888	DX-D 10G CASSETTE DETECTOR DRZ+4336R WITHOUT ENCL	-
2	RF+39889	DX-D 20G PORTABLE DETECTOR, RECEPTOR DRZ+ 4336R REV. A AND B - (REP.)	-
2	VA+39889	DX-D 20G PORTABLE DETECTOR DRZ+4336R WITH ENCL	-
3	RF+29886	DX-D 10C CASSETTE DETECTOR, RECEPTOR CSI 4336R REV. A AND B - (REP.)	-
3	VA+29886	DX-D 10C CASSETTE DETECTOR CSI 4336R WITHOUT ENCL.	-
4	RF+29887	DX-D 20C PORTABLE DETECTOR, RECEPTOR CSI 4336R REV. A AND B - (REP.)	-
4	VA+29887	DX-D 20C PORTABLE DETECTOR CSI 4336R WITH ENCL.	-
5	VA+13422	IO BOX (4336 V2)	-
6	VA+21638	CABLE, PIGTAIL 8M	-
7	VA+12792	CABLE, PIGTAIL 1M	-
8	VA+17463	CABLE, EXTENSION 7M	-
9	VA+32073	CABLE, PRIMARY, 4.5M HOT SWAPPABLE (V2) FOR PANEL - WITH ENCL.	-
10	VA+32075	CABLE, EXTENSION 3.5M	-
11	VA+17461	CABLE, EXTENSION 2M HOT SWAPPABLE	-
12	VA+660	CORD, MAINS, EUROPEAN	-
99	VA+15644	GRID HOLDER	-
99	VA+19751	OVERLAY DX-D 10 C	-
99	VA+19757	OVERLAY DX-D 20 G	-
99	VA+20552	OVERLAY DX-D 10 G	-
99	VA+20562	OVERLAY DX-D 20 C	-
99	VA+2517	AGFA 4336R (ENCLOSURE HOUSING)	-
99	VA+35291	CABLE, PRIMARY, 4.5M HOT SWAPPABLE (V2) FOR PANEL - WITHOUT ENCL.	-

□ = Assembly

* = Spare Part standardly is not kept in stock, extended delivery time is possible.

99 = Spare Part is not shown in Exploded Views.

MOBILE X-RAY UNIT


DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.


01



02



03 - 04



MOBILE X-RAY UNIT

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.

Item No.	Part No.	Description	Replacement Instruction ID
1	CM+6076350	INTEL GIGABIT CT DESKTOP CARD	-
2	CM+6213690	INTEL ETHERNET I210-T1 GBE NIC CARD (HP)	-
3	10+9999913390	GRID DX-D 20 R10:1 215L/INCH F40-72" (5UNHT)	-
4	10+9999913030	GRID DX-D 20, R6:1, 215 L/INCH, F 40-72 (5QP8M)	-

□ = Assembly

* = Spare Part standardly is not kept in stock, extended delivery time is possible.

99 = Spare Part is not shown in Exploded Views.

MOBILE X-RAY UNIT

DOCUMENT CONTROL NOTE:

The controlled version of this document is available from the Agfa Medimg Library. Any printed copy of this document is uncontrolled.